

# INNO Implant System™: Compatibility Table

Submerged(Bone Level)				Internal(Tissue Level)				External(Bone Level)			
COWELLMEDI				COWELLMEDI				COWELLMEDI			
Product	Diameter(Ø)	Connection		Product	Diameter(Ø)	Connection		Product	Diameter(Ø)	Connection	
INNO Submerged	3.1/3.3	Hex 2.1	Narrow	INNO Internal	3.5(P4.8)	Octa 3.1		INNO External	3.5(P4.2)	Hex 2.7	
	3.5/4.0	Hex 2.5	Regular & Short		4.0(P4.8)	Octa 3.1			4.0(P4.2)	Hex 2.7	
	4.5	Hex 2.5			4.5(P4.8)	Octa 3.1			4.5(P4.6)	Hex 2.7	
	5.0	Hex 2.5			5.0(P5.9)	Octa 3.1			5.0(P5.1)	Hex 3.4	
	6.0	Hex 2.5			6.0(P5.9)	Octa 3.1			6.0(P6.0)	Hex 3.4	
Astra				Zimmer				3i Biomet			
Product	Diameter(Ø)	Connection	Compatibility	Product	Diameter(Ø)	Connection	Compatibility	Product	Diameter(Ø)	Connection	Compatibility
OsseoSpeed	3.0	X-Small	N	SwissPlus	3.7(P3.8)	Octa 2.7	N	Osseotite	3.25(P3.5)	Hex 2.4	N
	3.5	Small	N		3.7(P4.8)	Octa 3.1	Y		4.0(P4.1)	Hex 2.7	Y
	4.0	Small	N		4.1(P4.8)	Octa 3.1	Y		5.0(P5.0)	Hex 2.7	Y
	4.5	Large	Y		4.8(P4.8)	Octa 3.1	Y		6.0(P5.0)	Hex 2.7	Y
	5.0	Large	Y								
Osstem				Straumann				Nobelbiocare			
Product	Diameter(Ø)	Connection	Compatibility	Product	Diameter(Ø)	Connection	Compatibility	Product	Diameter(Ø)	Connection	Compatibility
GS II, GS III TS II, TS III, TS IV	3.0/3.5	Hex 2.1	Y	Roxolid SLActive	3.3(P3.5/P4.8)	Octa 3.1	Y	Mk III/Mk IV/ NobelSpeedy Groovy	3.3(P3.5)	Hex 2.4	N
	4.0	Hex 2.5	Y		4.1(P4.8)	Octa 3.1	Y		3.75(P4.1)	Hex 2.7	Y
	4.5	Hex 2.5	Y		4.8(P4.8/P6.5)	Octa 3.1	Y		4.0(P4.1)	Hex 2.7	Y
	5.0	Hex 2.5	Y						5.0(P5.1)	Hex 3.4	Y
	5.5/6.0/6.0/7.0	Hex 2.5	Y						6.0(P5.1)	Hex 3.4	Y
Dentium				Megagen				Osstem			
Product	Diameter(Ø)	Connection	Compatibility	Product	Diameter(Ø)	Connection	Compatibility	Product	Diameter(Ø)	Connection	Compatibility
Implantium/Super Line I, II	3.7(3.6)	Hex 2.5	Y	AnyOne (Tissue level)	3.5(P4.8)	Octa 3.1	Y	US II, US III	3.5(P3.5)	Hex 2.4	N
	4.0	Hex 2.5	Y		4.0(P4.8)	Octa 3.1	Y		4.0(P4.1)	Hex 2.7	Y
	4.5	Hex 2.5	Y		4.5(P4.8)	Octa 3.1	Y		4.5(P4.1)	Hex 2.7	Y
	5.0	Hex 2.5	Y		5.0(P4.8)	Octa 3.1	Y		5.0(P5.0)	Hex 2.7	N
	6.0/7.0	Hex 2.5	Y						5.0(P5.1)	Hex 3.4	Y
DIO				Neobiotech				DIO			
Product	Diameter(Ø)	Connection	Compatibility	Product	Diameter(Ø)	Connection	Compatibility	Product	Diameter(Ø)	Connection	Compatibility
UF/Uf II/Steady	3.8	Hex 2.5	Y	IT, IT II, IT III	3.5(P4.8)	Octa 3.1	Y	FTN/FSN	3.3(P3.5)	Hex 2.4	N
	4.0	Hex 2.5	Y		4.0(P4.8)	Octa 3.1	Y		3.8(P4.1)	Hex 2.7	Y
	4.5	Hex 2.5	Y		4.5(P4.8)	Octa 3.1	Y		4.0(P4.1)/5.0(P5.0)	Hex 2.7	Y
	5.0	Hex 2.5	Y		5.0(P6.5)	Octa 3.1	N		5.0(P5.1)	Hex 3.4	Y
	5.5/6.0/6.5/7.0	Hex 2.5	Y		6.0(P6.5)	Octa 3.1	N		5.5/6.0(P5.1)	Hex 3.4	Y
Megagen				Dentis				Megagen			
Product	Diameter(Ø)	Connection	Compatibility	Product	Diameter(Ø)	Connection	Compatibility	Product	Diameter(Ø)	Connection	Compatibility
AnyOne(Normal thread, Deep thread)	3.5	Hex 2.5	Y	Oneq (i-clean)	3.7(P4.8)	Octa 3.1	Y	AnyOne	3.5(P3.5)	Hex 2.4	N
	4.0	Hex 2.5	Y		4.2(P4.8)	Octa 3.1	Y		4.0(P4.1)	Hex 2.7	Y
	4.5	Hex 2.5	Y		4.7(P4.8/P6.5)	Octa 3.1	Y		4.5(P4.5)	Hex 2.7	Y
	5.5/6.0/6.5	Hex 2.5	Y		5.2(P4.8/P6.5)	Octa 3.1	N		5.0(P5.0)	Hex 2.7	Y
	7.0/7.5/8.0	Hex 2.5	Y		6.0/7.0(P6.5)	Octa 3.1	N		6.0(P5.5)	Hex 2.7	Y
Neobiotech				SIS				Neobiotech			
Product	Diameter(Ø)	Connection	Compatibility	Product	Diameter(Ø)	Connection	Compatibility	Product	Diameter(Ø)	Connection	Compatibility
IS/IS II/IS III active	3.5	Hex 2.5	Y	Luna	3.5(P3.5)	Hex 2.4	N	EB/Eb II/ EB S-Wide	3.5(P3.5)	Hex 2.4	N
	4.0	Hex 2.5	Y		4.0	Hex 2.5	Y		4.0(P4.1)	Hex 2.7	Y
	4.5	Hex 2.5	Y		4.5	Hex 2.5	Y		4.5(P4.1)/5.0(P5.0)	Hex 2.7	Y
	5.0	Hex 2.5	Y		5.0	Hex 2.5	Y		5.0(P5.1)	Hex 3.4	Y
	5.5/6.0	Hex 2.5	Y		5.5/6.0	Hex 2.5	Y		5.5/6.0(P5.1)	Hex 3.4	Y
SIS				Dentis				SIS			
Product	Diameter(Ø)	Connection	Compatibility	Product	Diameter(Ø)	Connection	Compatibility	Product	Diameter(Ø)	Connection	Compatibility
Oneq (S-clean)	3.7	Hex 2.5	Y	Oneq (i-clean)	3.3/3.5(P3.5)	Hex 2.4	N	Sola	3.3/3.5(P3.5)	Hex 2.4	N
	4.2	Hex 2.5	Y		4.0	Hex 2.5	Y		4.0(P4.1)	Hex 2.7	Y
	4.7	Hex 2.5	Y		4.5	Hex 2.5	Y		4.5(P4.1)	Hex 2.7	Y
	5.2	Hex 2.5	Y		5.0	Hex 2.5	Y		5.0(P5.1)	Hex 3.4	Y
	6.0/7.0/8.0	Hex 2.5	Y		6.0	Hex 2.5	Y		6.0(P5.1)	Hex 3.4	Y
IBS				Dentis				Dentis			
Product	Diameter	Connection	Compatibility	Product	Diameter(Ø)	Connection	Compatibility	Product	Diameter(Ø)	Connection	Compatibility
Magic FC/ Magic FC II	4.0	Hex 2.5	Y	Oneq (S-clean)	3.5(P3.5)	Hex 2.4	N	e-clean	3.5(P3.5)	Hex 2.4	N
	4.5	Hex 2.5	Y		4.1	Hex 2.5	Y		4.1(P4.1)	Hex 2.7	Y
	5.0	Hex 2.5	Y		4.7	Hex 2.5	Y		5.1(P5.1)	Hex 3.4	Y
	5.5	Hex 2.5	Y		5.2	Hex 2.5	Y				
	6.0/6.5	Hex 2.5	Y		6.0/7.0/8.0	Hex 2.5	Y				
Snucone				AF/AF+B/AF+I							
Product	Diameter	Connection	Compatibility	Product	Diameter	Connection	Compatibility				
AF/AF+B/AF+I	3.5	Hex 2.5	Y	AF/AF+B/AF+I	3.5	Hex 2.5	Y				
	4.0	Hex 2.5	Y		4.0	Hex 2.5	Y				
	4.3	Hex 2.5	Y		4.3	Hex 2.5	Y				
	4.8	Hex 2.5	Y		4.8	Hex 2.5	Y				
	5.3/5.8	Hex 2.5	Y		5.3/5.8	Hex 2.5	Y				